

Anti-Aldolase Rabbit antibody



Catalog Number: 167214

Product name

Rabbit IgG

Anti-Aldolase Rabbit antibody

Applications

WB, IHC-P, ICC/IF, FC

Specificity

Human, Mouse, Rat

Dilutions

WB: 1/1000-1/5000

IHC: 1/20

ICC/IF: 1/20

FC: 1/50

Antibody description

Aldolase Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Aldolase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Validations



Western blot detection of Aldolase in K562, Rat Brain, C6, 3T3, HeLa cell lysates using Aldolase Rabbit pAb(1:1000 diluted). Predicted band size: 39kDa. Observed band size: 39kDa.