

UNC45B antibody



Catalog Number: 116614

Product name

UNC45B antibody

Clonality

Polyclonal

Immunogen

[Human UNC45B Recombinant protein \(GST tag & His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

UNC45B Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 95kd

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of UNC45B recombinant protein (Accession Number: NM_001267052). Purification method: Antigen affinity purified.


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116614(UNC45B antibody) at dilution of 1:500