

Perforin antibody



Catalog Number: 113714

Product name

Perforin antibody

Immunogen

[Human Perforin Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

Perforin Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse liver tissue. Observed molecular weight by Western-blot: 70-75 kDa

Preparation

This antibody was obtained by immunization of Perforin recombinant protein (Accession Number: NM_001083116). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113714(PRF1 antibody) at dilution of 1:500