

# STK25 antibody



Catalog Number: 115714

## Product name

STK25 antibody

## Immunogen

[Human STK25 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

STK25 Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells, mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of STK25 recombinant protein (Accession Number: NM\_006374). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:115714(STK25 Antibody) at dilution of 1:600