

Anti-Ngfr antibody



Catalog Number: 175024

Product name

Anti-Ngfr antibody

Specificity

Mouse, Rat

Antibody description

Rabbit polyclonal antibody to Ngfr

Preparation

This antigen of this antibody was peptide

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC, IHC-P, FC

Dilutions

WB: 1:2,000~1:5,000

Validations




Fig1: Western blot analysis on D3 (A) and PC12 (B) using anti-CD271 polyclonal antibody.