

Anti-Aryl hydrocarbon Receptor Rabbit antibody



Catalog Number: 166824

Product name WB, ICC/IF

Anti-Aryl hydrocarbon Receptor Rabbit antibody **Dilutions**

Specificity WB: 1/1000

Human ICC/IF: 1/20

Antibody description

Aryl hydrocarbon Receptor Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Aryl hydrocarbon Receptor

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



Immunocytochemistry of Aryl hydrocarbon Receptor(green) in HeLa cells using Aryl hydrocarbon Receptor Rabbit pAb at dilution 1/50,



and DAPI(blue) 17-
Western blot detection of Aryl hydrocarbon Receptor in HeLa cell lysates using Aryl hydrocarbon Receptor Rabbit pAb(1:1000 diluted).Predicted band size:96kDa.Observed band size:100kDa.