

# HOPX antibody



Catalog Number: 111524

## Product name

HOPX antibody

## Immunogen

[Human HOPX Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HOPX Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human placenta tissue, rat lung tissue. Positive IP detected in Mouse lung. Observed molecular weight by Western-blot: 8-12kd

## Preparation

This antibody was obtained by immunization of HOPX recombinant protein (Accession Number: NM\_139212). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



IP result of HOPX antibody (Catalog No:111524 for IP and Detection) with mouse lung tissue



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111524(HOPX Antibody) at dilution of 1:600