

# Anti-ULK1 Rabbit antibody

Catalog Number: 167524

## Product name

Anti-ULK1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

ULK1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human ULK1

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/500-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200


## Validations



Western blot analysis of ULK1 expression in (1) HEK293 cell lysate; (2) PC12 cell lysate.



Immunohistochemical analysis of paraffin-embedded human heart, using ULK1 Antibody.



Immunofluorescent analysis of 293 cells, using ULK1 Antibody .