

VBP1 antibody



Catalog Number: 116724

Product name

VBP1 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VBP1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive IF detected in HepG2 cells. Positive WB detected in rat brain tissue, C6 cells, HepG2 cells, K-562 cells, mouse liver tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 23-25 kDa

Preparation

This antibody was obtained by immunization of VBP1 recombinant protein (Accession Number: NM_003372). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IF, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116724(VBP1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116724(VBP1 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using VBP1 antibody Catalog No:116724 at 1:50 dilution

VBP1 antibody



Catalog Number: 116724

and Rhodamine-labeled goat anti-rabbit IgG (red).
Blue pseudocolor = DAPI (fluorescent DNA dye).

rat brain tissue were subjected to SDS PAGE
followed by western blot with Catalog
No:116724(VBP1 Antibody) at dilution of 1:1000

