

# ACAP1 antibody



Catalog Number: 107724

## Product name

ACAP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ACAP1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human lung tissue. Observed molecular weight by Western-blot: 85 kDa

## Preparation

This antibody was obtained by immunization of ACAP1 recombinant protein (Accession Number: NM\_014716). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107724(ACAP1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung using Catalog No:107724(ACAP1 antibody) at dilution of 1:100 (under 10x lens)