

MORF4L2 Polyclonal Antibody



Catalog Number: 165434

Product name

MORF4L2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MORF4L2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MORF4L2 (NP_001135890.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:100

Validations



Immunofluorescence - MORF4L2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using MORF4L2 antibody .