

EMR1 Polyclonal Antibody



Catalog Number: 161134

Product name

EMR1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to EMR1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-290 of human EMR1 (NP_001243182.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - EMR1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using EMR1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - EMR1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using EMR1 antibody at dilution of



1:200 (40x lens).

EMR1 Polyclonal Antibody



Catalog Number: 161134

Immunohistochemistry - EMR1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spleen using EMR1 antibody at dilution of 1:200 (40x lens).