

Product name

Anti-NOL8 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to NOL8

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human nol8. aa range:1118-1167

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300


Validations

Fig1:: Western Blot analysis of 3T3 cells using Nop132 Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).