

Anti-TXNL1 Rabbit antibody



Catalog Number: 167234

Product name

Anti-TXNL1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TXNL1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human TXNL1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/5000

IHC: 1/20


ICC/IF: 1/20

FC: 1/20


Validations



Immunohistochemistry of TXNL1 in paraffin-embedded Human breast cancer tissue using TXNL1 Rabbit pAb at dilution 1/20



Western blot detection of TXNL1 in K562,Rat Brain,C6,3T3,Hela cell lysates using TXNL1 Rabbit pAb(1:1000 diluted).Predicted band size:32kDa.Observed band size:32kDa.



Immunocytochemistry of TXNL1(green) in HeLa cells using TXNL1 Rabbit pAb at dilution 1/50, and DAPI(blue)