

TPPII antibody



Catalog Number: 107634

Product name

TPPII antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

[Human TPP2 Recombinant protein \(His tag\)](#)


IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human; other species not tested.

Validations



Antibody description

TPPII Mouse Monoclonal antibody. Positive FC detected in HeLa cells. Positive IHC detected in human testis tissue. Positive IF detected in HeLa cells, MCF-7 cells. Positive WB detected in HeLa cells, HepG2 cells, Raji cells. Observed molecular weight by Western-blot: 135 kDa

Preparation

This antibody was obtained by immunization of TPPII recombinant protein (Accession Number: NM_003291). Purification method: Protein A purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107634(TPP2 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal


Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC, IF, FC

Dilutions




Immunohistochemical of paraffin-embedded human testis using Catalog No:107634(TPP2 antibody) at dilution of 1:50 (under 10x lens)


TPPII antibody




Catalog Number: 107634



Immunohistochemical of paraffin-embedded human testis using Catalog No:107634(TPP2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using TPP2 antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).



1×10^6 HeLa cells were stained with 0.2ug TPP2 antibody (Catalog No:107634, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.