

HAL antibody



Catalog Number: 111334

Product name

HAL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HAL Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, rat liver tissue. Positive IHC detected in human skin tissue. Observed molecular weight by Western-blot: 73 kDa

Preparation

This antibody was obtained by immunization of HAL recombinant protein (Accession Number: NM_002108). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111334(HAL Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:111334(HAL Antibody) at dilution of 1:50 (under 40x lens)