

# PFDN2 antibody



Catalog Number: 113734

## Product name

PFDN2 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

## Specificity

Human; other species not tested.

## Antibody description

PFDN2 Rabbit Polyclonal antibody. Positive WB detected in L02 cells. Positive IHC detected in human lung cancer tissue. Positive IF detected in L02 cells. Observed molecular weight by Western-blot: 17-20 kDa

## Preparation

This antibody was obtained by immunization of PFDN2 recombinant protein (Accession Number: NM\_012394). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000


## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:113734(PFDN2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113734(PFDN2 antibody) at dilution of 1:100 (under 10x lens)



Immunofluorescent analysis of L02 cells using

# PFDN2 antibody



Catalog Number: 113734

---

Catalog No:113734(PFDN2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)