

CD14 antibody



Catalog Number: 107044

Product name

CD14 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD14 Mouse Monoclonal antibody. Positive WB detected in RAW264.7 cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of CD14 recombinant protein (Accession Number: NM_001174104). Purification method: .

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107044(CD14 antibody) at dilution of 1:500