

MRPS17 antibody



Catalog Number: 112844

Product name

MRPS17 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

[Human MRPS17 Recombinant protein \(His tag\)](#)

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

MRPS17 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Positive IHC detected in human kidney tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 15 kDa




Preparation

This antibody was obtained by immunization of MRPS17 recombinant protein (Accession Number: NM_015969). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112844(MRPS17 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human kidney using Catalog No:112844(MRPS17 antibody) at dilution of 1:50 (under 10x lens)


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions




MRPS17 antibody



Catalog Number: 112844

Immunohistochemical of paraffin-embedded human kidney using Catalog No:112844(MRPS17 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:112844(MRPS17 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)