

NFE2 antibody



Catalog Number: 113144

Product name

NFE2 antibody

Immunogen

[Human NFE2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

NFE2 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive IF detected in HepG2 cells. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, U-937 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of NFE2 recombinant protein (Accession Number: NM_001136023). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113144(NFE2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113144(NFE2 antibody) at dilution of 1:50 (under 10x lens)

NFE2 antibody

Catalog Number: 113144



Immunofluorescent analysis of HepG2 cells, using NFE2 antibody Catalog No:113144 at 1:100 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent

DNA dye).

IP Result of anti-NFE2 (IP:Catalog No:113144, 4ug; Detection:Catalog No:113144 1:500) with K-562 cells lysate 3440ug.

