

MOK Polyclonal Antibody



Catalog Number: 163544

Product name

MOK Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MOK

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human MOK (NP_055041.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - MOK Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using MOK antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - MOK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using MOK antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - MOK Polyclonal Antibody

MOK Polyclonal Antibody



Catalog Number: 163544

Immunohistochemistry of paraffin-embedded rat heart using MOK antibody at dilution of 1:100




(40x lens).

Immunohistochemistry - MOK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded

human liver using MOK antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - MOK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using MOK antibody at dilution of 1:100 (40x lens).