

Anti-beta Dystroglycan Rabbit antibody



Catalog Number: 167944

Product name

Anti-beta Dystroglycan Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

beta Dystroglycan Rabbit monoclonal antibody


Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic peptide of human beta Dystroglycan

Validations



Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of beta Dystroglycan in MCF7,HepG2,Lncap cell lysates using beta Dystroglycan Rabbit mAb(1:1000 diluted).Predicted band size:27kDa.Observed band size:42-44kDa.