

SPRR2D antibody



Catalog Number: 115644

Product name

SPRR2D antibody

Specificity

Human; other species not tested.

Antibody description

SPRR2D Rabbit Polyclonal antibody. Positive IHC detected in human oesophagus tissue, human skin tissue.

Preparation

This antibody was obtained by immunization of SPRR2D recombinant protein (Accession Number: XM_017002175). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:115644 (SPRR2D Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:115644 (SPRR2D Antibody) at dilution of 1:50 (under 40x lens)