

URM1 antibody



Catalog Number: 116644

Product name

URM1 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

IHC: 1:20-1:200

[Human URM1 Recombinant protein \(GST tag & His tag\)](#)


Validations

Specificity

human, mouse, rat; other species not tested.

Antibody description

URM1 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, A549 cells, HepG2 cells, mouse liver tissue. Positive IHC detected in human hepat°Cirrhosis. Observed molecular weight by Western-blot: 11kd




Preparation

This antibody was obtained by immunization of URM1 recombinant protein (Accession Number: NM_030914). Purification method: Antigen affinity purified.

human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116644(URM1 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116644(URM1 antibody) at dilution of 1:100 (under 10x lens)

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions



URM1 antibody



Catalog Number: 116644

Immunohistochemical of paraffin-embedded
human hepatocirrhosis using Catalog

No:116644(URM1 antibody) at dilution of 1:100
(under 40x lens)