

PVRL1 antibody



Catalog Number: 114344

Product name

PVRL1 antibody

IHC: 1:50-1:500

IF: 1:50-1:500

Specificity

Human, Mouse; other species not tested.

Antibody description

PVRL1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Positive IF detected in SH-SY5Y cells. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of PVRL1 recombinant protein (Accession Number: NM_002855). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:114344(PVRL1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114344(PVRL1 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).




PVRL1 antibody



Catalog Number: 114344

Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114344(PVRL1 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunofluorescent analysis of (-20oc Ethanol) fixed SH-SY5Y cells using Catalog No:114344(PVRL1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)