

SAMD9 antibody



Catalog Number: 115054

Product name

SAMD9 antibody

Clonality

Polyclonal

Immunogen

[Human SAMD9 Recombinant protein \(GST tag & His tag\)](#)

Ig Type

Rabbit IgG

Specificity

human; other species not tested.

Applications

ELISA, WB

Antibody description

SAMD9 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 148 kDa


Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of SAMD9 recombinant protein (Accession Number: NM_017654). Purification method: Antigen Affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115054(SAMD9 antibody) at dilution of 1:1000