

Product name

NXN antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Validations**Antibody description**

NXN Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human heart tissue. Positive IP detected in HeLa cells. Positive IF detected in HepG2 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 48 kDa

**Preparation**

This antibody was obtained by immunization of NXN recombinant protein (Accession Number: NM_001205319). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113454(NXN antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF, IP


Dilutions

Recommended Dilution:

WB: 1:200-1:2000




Immunohistochemical of paraffin-embedded human testis using Catalog No:113454(NXN antibody) at dilution of 1:50 (under 10x lens)



NXN antibody


Catalog Number: 113454

Immunohistochemical of paraffin-embedded human testis using Catalog No:113454(NXN antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using

NXN antibody Catalog No:113454 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



IP Result of anti-NXN (IP:Catalog No:113454, 4ug; Detection:Catalog No:113454 1:500) with HeLa cells lysate 3200ug.