

PKD2L1 antibody



Catalog Number: 113854

Product name

Polyclonal

PKD2L1 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

PKD2L1 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 92kd

Dilutions

Recommended Dilution:

Preparation


WB: 1:200-1:2000

This antibody was obtained by immunization of PKD2L1 recombinant protein (Accession Number: NM_016112). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113854(PKD2L1 antibody) at dilution of 1:500