

# Anti-APOBEC3D antibody



Catalog Number: 175854

## Product name

Polyclonal

Anti-APOBEC3D antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Monkey

## Applications

## Antibody description

WB, IHC-P

Rabbit polyclonal antibody to APOBEC3D

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human apobec3d/f. aa range:232-281

IHC-P: 1:100-1:300

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

## Validations




Fig1.; Western Blot analysis of various cells using APOBEC3D/F Polyclonal Antibody