

Anti-GNAQ Rabbit antibody



Catalog Number: 168854

Product name

Anti-GNAQ Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

GNAQ Rabbit monoclonal antibody

Dilutions

1:1000

Preparation

Antigen: A synthetic peptide of human GNAQ

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of GNAQ in Rat Brain, C6, 3T3 cell lysates using GNAQ Rabbit mAb(1:1000 diluted). Predicted band size: 42kDa. Observed band size: 42kDa.