

EFTUD2 Polyclonal Antibody



Catalog Number: 165154

Product name

EFTUD2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to EFTUD2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human EFTUD2 (NP_004238.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF, IP

Dilutions


WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:100


IP 1:20 - 1:50

Validations




Western blot - EFTUD2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using EFTUD2 antibody at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 20s.



Immunohistochemistry - EFTUD2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using EFTUD2 Antibody at dilution of



1:100 (40x lens).


EFTUD2 Polyclonal Antibody



Catalog Number: 165154

Immunohistochemistry - EFTUD2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat liver using EFTUD2 Antibody at dilution of 1:100



(40x lens).


Immunohistochemistry - EFTUD2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using EFTUD2 Antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - EFTUD2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using EFTUD2 Antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - EFTUD2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using EFTUD2 Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - EFTUD2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using EFTUD2 Antibody at dilution



Immunofluorescence - EFTUD2 Polyclonal Antibody

EFTUD2 Polyclonal Antibody



Catalog Number: 165154

Immunofluorescence analysis of U2OS cells using EFTUD2 antibody .