

# ZNF277 antibody



Catalog Number: 117154

## Product name

ZNF277 antibody

## Immunogen

[Human ZNF277 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

ZNF277 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in PC-3 cells, Jurkat cells. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of ZNF277 recombinant protein (Accession Number: BC020626). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:117154(ZNF277 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:117154(ZNF277 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG