

ARMCX2 antibody



Catalog Number: 108154

Product name

ARMCX2 antibody

Clonality

Polyclonal

Immunogen

[Human ARM CX2 Recombinant protein \(GST tag & His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

ARMCX2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse brain tissue, mouse liver tissue, mouse testis tissue. Observed molecular weight by Western-blot: 66 kDa

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of ARM CX2 recombinant protein (Accession Number: XM_047442689). Purification method: Antigen affinity purified.

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108154(ARM CX2 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT