

FBXO41 antibody



Catalog Number: 110554

Product name

FBXO41 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FBXO41 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human brain tissue. Positive WB detected in mouse brain tissue, rat brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of FBXO41 recombinant protein (Accession Number: BC112202). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110554(FBXO41 Antibody) at dilution of



IP Result of anti-FBXO41 (IP:Catalog No:110554, 4ug; Detection:Catalog No:110554 1:1000) with mouse brain tissue lysate 6000ug.




Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog

FBXO41 antibody



Catalog Number: 110554

No:110554(FBXO41 Antibody) at dilution of 1:50
(under 40x lens)



Immunohistochemistry of paraffin-embedded
human heart tissue slide using Catalog
No:110554(FBXO41 Antibody) at dilution of 1:50
(under 40x lens)