

PLEK2 antibody



Catalog Number: 113954

Product name

PLEK2 antibody

Immunogen

[Human PLEK2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLEK2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue, mouse spleen tissue. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of PLEK2 recombinant protein (Accession Number: NM_016445). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113954(PLEK2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113954(PLEK2 antibody) at dilution of 1:50 (under 10x lens)