

Anti-TSFM Rabbit antibody



Catalog Number: 167954

Product name

Rabbit IgG

Anti-TSFM Rabbit antibody

Applications

Specificity

WB, ICC/IF

Human

Dilutions

Antibody description

WB: 1/1000

TSFM Rabbit monoclonal antibody

ICC/IF: 1/20

Preparation

Antigen: A synthetic peptide of human TSFM

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of TSFM in HeLa,A549,HL-60,U2OS cell lysates using TSFM Rabbit mAb(1:1000 diluted).Predicted band size:35kDa.Observed band size:35kDa.