

# Anti-ERCC1 Rabbit antibody



Catalog Number: 168654

## Product name

WB, ICC/IF, FC

Anti-ERCC1 Rabbit antibody

## Dilutions

## Specificity

WB: 1/1000-1/5000

Human, Mouse

ICC/IF: 1/20-1/100

## Antibody description

FC: 1/50

ERCC1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human ERCC1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


Rabbit IgG

## Applications

## Validations



Immunocytochemistry of ERCC1(green) in HeLa cells using ERCC1 Rabbit mAb at dilution 1/200,



Western blot detection of ERCC1 in K562,3T3,HeLa cell lysates using ERCC1 Rabbit mAb(1:1000 diluted).Predicted band size:33kDa.Observed band size:39kDa.