

Human DTD1/C20orf88/DUEB/HARS2 (His tag) recombinant protein



Catalog Number: 519754

General Information

Gene Name Synonym

D-tyrosyl-tRNA(Tyr) deacylase 1, DNA-unwinding element-binding protein B, DUE-B, Histidyl-tRNA synthase-related, DTD1, C20orf88, DUEB, HARS2

Protein Construction

A DNA sequence encoding the Human DTD1/C20orf88/DUEB/HARS2 Met1-Pro209 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human DTD1/C20orf88/DUEB/HARS2 (His tag)
recombinant protein