

POSH antibody



Catalog Number: 114064

Product name

POSH antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

POSH Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse brain tissue, DU 145 cells, HEK-293 cells, human brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 93-100 kDa

Preparation

This antibody was obtained by immunization of POSH recombinant protein (Accession Number: BC053671). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114064(SH3RF1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114064(SH3RF1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

POSH antibody



Catalog Number: 114064

human gliomas using Catalog No:114064(SH3RF1 antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-SH3RF1 (IP:Catalog No:114064, 4ug; Detection:Catalog No:114064 1:1000) with mouse brain tissue lysate 4000ug.

