

APRT Polyclonal Antibody



Catalog Number: 161864

Product name

APRT Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to APRT

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human APRT (NP_000476.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - APRT Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using APRT antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry - APRT Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using APRT Antibody at dilution of 1:100 (40x lens).



APRT Polyclonal Antibody



Catalog Number: 161864

Immunohistochemistry - APRT Polyclonal Antibody mouse lung using APRT Antibody at dilution of
Immunohistochemistry of paraffin-embedded 1:100 (40x lens).