

Anti-Cyclin B1 Rabbit antibody



Catalog Number: 167164

Product name

WB, ICC/IF, FC

Anti-Cyclin B1 Rabbit antibody

Dilutions

Specificity

WB: 1/1000

Human, Mouse

ICC/IF: 1/100

Antibody description

FC: 1/50

Cyclin B1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of mouse Cyclin B1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



Immunocytochemistry of Cyclin B1(green) in HeLa cells using Cyclin B1 Rabbit pAb at dilution 1/50,



and DAPI(blue)

Western blot detection of Cyclin B1 in Jurkat cell lysates using Cyclin B1 Rabbit pAb(1:1000 diluted). Predicted band size: 48KDa. Observed band size: 55KDa.