

C10orf46 antibody



Catalog Number: 108564

Product name

Polyclonal

C10orf46 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

C10orf46 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, L02 cells. Observed molecular weight by Western-blot: 41 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of C10orf46 recombinant protein (Accession Number: BC038294). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108564(C10orf46 Antibody) at dilution of 1:600