

# TAB1 antibody



Catalog Number: 115964

## Product name

TAB1 antibody

## Immunogen

[Human MAP3K7IP1 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TAB1 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of TAB1 recombinant protein (Accession Number: NM\_006116). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:115964(MAP3K7IP1 antibody) at dilution of



IP Result of anti-MAP3K7IP1 (IP:Catalog No:115964, 4ug; Detection:Catalog No:115964 1:600) with K-562 cells lysate 3200ug.