

Anti-NRF1 Rabbit antibody



Catalog Number: 168964

Product name 1:2000

Anti-NRF1 Rabbit antibody 1:20

Specificity 1:20

Human, Mouse, Rat 1:20

Antibody description 1:20

NRF1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NRF1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of NRF1 in C6,3T3,Hela cell lysates using NRF1 Rabbit mAb(1:1000 diluted).Predicted band size:54kDa.Observed band size:68kDa.



Immunocytochemistry of NRF1(green) in HeLa cells using NRF1 Rabbit mAb at dilution 1/200, and DAPI(blue)