

PRL Polyclonal Antibody



Catalog Number: 163074

Product name

PRL Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to PRL

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 29-227 of human PRL (NP_001157030.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:10 - 1:100

Validations



Western blot - PRL Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PRL antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - PRL Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using PRL antibody .