

TBCEL antibody



Catalog Number: 115874

Product name

TBCEL antibody

Immunogen

[Human TBCEL Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TBCEL Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, PC-3 cells. Observed molecular weight by Western-blot: 40-48 kDa

Preparation

This antibody was obtained by immunization of TBCEL recombinant protein (Accession Number: NM_001363644). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115874(TBCEL Antibody) at dilution of 1:3000