

Zyxin antibody



Catalog Number: 117274

Product name

Zyxin antibody

Immunogen

[Human ZYX Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

Zyxin Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF7 cells. Observed molecular weight by Western-blot: 78kd

Preparation

This antibody was obtained by immunization of Zyxin recombinant protein (Accession Number: NM_001010972). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117274(ZYX antibody) at dilution of 1:1000