

MAFF antibody



Catalog Number: 112374

Product name

MAFF antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MAFF Rabbit Polyclonal antibody. Positive WB detected in human lung tissue. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 18-23 kDa

Preparation

This antibody was obtained by immunization of MAFF recombinant protein (Accession Number: NM_012323). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112374(MAFF antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:112374(MAFF antibody) at dilution of 1:100 (under 10x lens)