

CBX8 Polyclonal Antibody



Catalog Number: 164774

Product name

CBX8 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CBX8

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 80-280 of human CBX8 (NP_065700.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:10 - 1:100

Validations




Western blot - CBX8 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CBX8 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry - CBX8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CBX8 Antibody at dilution of



1:200 (40x lens).


CBX8 Polyclonal Antibody



Catalog Number: 164774

Immunohistochemistry - CBX8 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using CBX8 Antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CBX8 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat pancreas using CBX8 Antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CBX8 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human colon using CBX8 Antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CBX8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse testis using CBX8 Antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CBX8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CBX8 Antibody at dilution of 1:100 (40x lens).