

Anti-Caspase-3 Rabbit antibody



Catalog Number: 167084

Product name

WB: 1/1000

Anti-Caspase-3 Rabbit antibody

IHC: 1/50

Specificity

IP: 1/20

Human, Mouse, Rat

Antibody description

Caspase-3 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Caspase-3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, IP

Dilutions

Validations



Western blot detection of Caspase-3 in K562,Rat Brain,C6,3T3,Hela cell lysates using Caspase-3 Rabbit pAb(1:1000 diluted).Predicted band size:32kDa.Observed band size:32kDa.



Immunohistochemistry of Caspase-3 in paraffin-embedded Human tonsil using Caspase-3 Rabbit pAb at dilution 1/100